



Sterile, Hydrated & Ready To Use.

**Every second counts.
Save time in the O.R.**

KeraLink's sterile scleral patch grafts arrive hydrated in saline, ready to use for your glaucoma drainage device procedures. Store these sterile patch grafts at room temperature and have them available for those unscheduled cases or if needed for tectonic support, revisions of trabeculectomy flaps or exposed tube shunts.

Sterile

- Gamma-irradiated
- No rehydration required
- No rinsing needed

Readily Available

- 2-Year shelf life
- Store at room temperature
- Cut to 1x1cm or 7x7mm patches

Innovative Tissue Processing

KeraLink's innovative process inactivates the antigenicity of cells without exposing patient to toxic materials. All grafts are terminally sterilized with precision gamma irradiation to a sterility assurance level of 10⁻⁶ and stored in balanced saline solution.

CATALOG NO.	DESCRIPTION	DIMENSIONS
S0500SI-11	Square Scleral Patch in Saline	1.0cm x 1.0cm
S0500SI-77	Square Scleral Patch in Saline	7.0mm x 7.0mm

Advancing Corneal Transplantation

For more than five decades, KeraLink International (formerly Tissue Banks International) has been a leader in corneal transplantation. Our singular mission is to reverse corneal blindness globally through transplantation, innovation, and partnerships. Much of our work is technical, but underlying all of our efforts is a deep respect for the remarkable selfless gift of tissue donation and a steadfast commitment to honor this gift.

KeraLink.org | 800.858.2020

5520 Research Park Drive
Suite 400
Baltimore, MD 21228
Tissue@KeraLink.org

Doc. No. OI-020 Effective Date: March 15, 2017

*Exclusively processed
and distributed by:*



KeraLink
INTERNATIONAL

FORMERLY
TBI/TISSUE BANKS INTERNATIONAL