

# Ocular Tissue Services



DMEK

Pre-stripped, pre-punched with I-II orientation mark.



VisionGraft®

Sterile half-moon cornea patch graft.



Sclera

Hydrated, sterile sclera patch graft.

## Advancing Corneal Transplantation

For more than five decades, KeraLink International (formerly Tissue Banks International) has been a leader in corneal transplantation. Today, KeraLink is made up of the largest network of non-profit, accredited eye banks nationwide. Our singular mission is to reverse corneal blindness globally through transplantation, innovation, and partnerships. Much of our work is technical, but within all of our efforts is a deep respect for the remarkable, selfless gift of tissue donation and a steadfast commitment to honor this gift.

## Live Cell Corneas


Evaluated and released for transplantation and stored in Opitsol-GS.

## Full Thickness/Pre-cut/Prepared Corneas

CATALOG NO.	DESCRIPTION	PREPARATION
K001-PK	Full thickness cornea	For PKP, procedures where graft is prepared by surgeon, or KLAL (limbal cells)
K002-DS	DSAEK	Custom prepared to thickness requested, including ultra-thin
K005-DM	DMEK	Pre-stripped w/ 1-2 mm center attachment. Pre-punched diameter (custom sized graft: 7.5 mm - 8.5 mm)
K003-AK	ALK/DALK	Full thickness or pre-cut customized thickness to surgeon's specifications

## VisionGraft®

Clear sterile cornea - gamma irradiated, stored in albumin at room temperature, ready to use to promote ease of use.



## For Glaucoma Procedures

CATALOG NO.	DESCRIPTION	DIMENSIONS
C0301AL-90	Split Thickness Half Moon	4.5 mm x 9.0 mm
C0302AL-90	Split Thickness Third Moon	4.5 mm x 6.75 mm (7.8 mm diag.)

## For Cornea Procedures

CATALOG NO.	DESCRIPTION	DIMENSIONS
C0100AL	Full Thickness Whole Moon	9.00 mm
C0101AL	Full Thickness Whole Moon with Scleral Rim	16.00 mm
C0300AL-90	Split Thickness Whole Moon	9.00 mm

## For Keratoprosthesis

CATALOG NO.	DESCRIPTION	DIMENSIONS
C0400AL-85	KPro Ring	8.5 mm, 3.0 mm center
C0400AL-90	KPro Ring	9.0 mm, 3.0 mm center

## Sterile Sclera

Gamma irradiated, stored in saline at room temperature, arrives hydrated and ready to use.

CATALOG NO.	DESCRIPTION	DIMENSIONS
S0500SI-11	Square Sclera Patch in Saline	1.0 cm x 1.0 cm
S0500SI-77	Square Sclera Patch in Saline	7.0 mm x 7.0 mm

**KeraLink.org | 800.858.2020**

5520 Research Park Drive  
Suite 400  
Baltimore, MD 21228  
Tissue@KeraLink.org

Doc. No. OA-016 Effective Date: September 15, 2017



**KeraLink™**  
INTERNATIONAL

FORMERLY  
**TBI/TISSUE BANKS INTERNATIONAL**